

**Catalog No.:** G915

**Product Name:** **Direct RT-qPCR Lysis Kit**

**Size:** 50 Preps

**Description:** **Direct RT-qPCR Lysis Kit** offers a simple and quick method to prepare template directly from  $10^5$  cultured cells. The kit includes reagents for cell lysis and gDNA removal at the same time. Without further RNA extraction and purification, the prepared template in the lysate can be used directly for reverse transcription and qPCR (two step RT-PCR); or directly for qRT-PCR (one step RT-PCR).

**Kit Contents:**

Components	Volume
Lysis Solution	2x 1.25ml
Stop Solution	300 $\mu$ l
Protease	50 $\mu$ l
Protease Inhibitor	50 $\mu$ l

**Storage:** Store all components at  $-20^{\circ}\text{C}$  in a non-defrosting freezer.

**Protocol:**

1. Thaw and leave the Lysis Solution and Stop Solution at room temperature and mix gently but thoroughly.
2. Prepare  $10^5$  cells per reaction for lysis.
  - **For adherent cells grown in 96- or 384-well plates:** Aspirate the culture medium from the wells and rinse with 50  $\mu$ l of chilled 1X PBS. Remove the PBS.
  - **For suspension cells:** Pellet up to  $10^5$  cells in a centrifuge tube. Remove the culture medium and rinse with 50  $\mu$ l of chilled 1X PBS. Remove the PBS.
3. **Prepare the Lysis Mix:** use a new centrifuge tube, mix 50  $\mu$ l Lysis Solution with 1  $\mu$ l Protease for each sample. For example, if you have 5 samples, you need to mix 250  $\mu$ l Lysis Solution with 5  $\mu$ l Protease.
4. **Prepare the Stop Mix:** use a centrifuge tube, mix 5  $\mu$ l Stop Solution with 1  $\mu$ l Protease Inhibitor for each sample; e.g., if you have 5 samples, you need to mix 25  $\mu$ l Stop Solution with 5  $\mu$ l Protease Inhibitor.
5. To lyse the cells, add 50  $\mu$ l of the **Lysis Mix** from Step 3 to the cells prepared from Step 2. Mix the cells well by pipetting 35  $\mu$ l of the mixture up and down for five times. Avoid creating bubbles.
6. Incubate for 10 minutes at  $37^{\circ}\text{C}$ .
7. To terminate the lysis, add 5  $\mu$ l of the **Stop Mix** prepared at Step 4; and mix the lysate by pipetting 35  $\mu$ l of the mixture up and down 5 times.
8. Incubate the final mixture for 5 minutes at room temperature. The lysates containing RNA templates is now ready for various downstream applications.

**Note:**

1. To minimize RNA degradation, keep the cells in PBS on ice before starting the cell lysis procedure.
2. Lysis Solution and Stop Solution must be at room temperature and mixed well before adding to the cells.
3. To protect the RNA templates in solution, we recommend using the lysates immediately within one hour. However, the lysates can be stored at  $-80^{\circ}\text{C}$  for short term storage with one freeze-thaw cycle.
4. The prepared lysate should be used less than 10% of the final volume for any downstream applications.