



**Catalog No.:** C108

**Product Name:** Yeast Extract

**Size:** 500g

**CAS#:** 8013-01-2

**Description:** High quality Yeast Extract. This Yeast Extract contains a mixture of vitamins, amino acids, peptides and carbohydrates. It is widely used for general bacteriological research work with a variety of microorganisms, such as for making various microbiological media and as an additive for other culture media.

**Specifications:** Typical representative sample of different lots of yeast extract.

Amino nitrogen:	5.4%
Total nitrogen:	10.3%
Appearance:	yellowish/off-white powder
Loss on drying:	2.8%

**Storage:** Room temperature or below, under dry conditions.

*Note: This Product Is For Research Use Only.*